

201.4 - Ion-Selective Electrode Calibration (powder form)

These SRMs are certified for the calibration of ion-selective electrodes and have conventional ionic activities based on the Stokes-Robinson hydration theory for ionic strengths greater than 0.1mol/L.

Technical Contact: kenneth.pratt@nist.gov

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

| SRM | 2201 | 2202 | 2203 |
|---------------------|---------------------------------|---|---|
| Description | Sodium Chloride (Ion-Selective) | Potassium Chloride (Ion-Selective Electrodes) | Potassium Fluoride (Ion-Selective Electrodes) |
| Unit of Issue | (125 g) | (160 g) | (125 g) |
| | | | |
| Certified Component | pNa, pCl | pK, pCl | pF |