203.1 - Combustion Calorimetry (powder form)

These SRMs are for use as standards for calibration of combustion bomb calorimeters used in checking the performance of apparatus and analytical procedures.

Technical Contact: donald.archer@nist.gov

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM	39 j	1656	1657	2151	2152
Description	Benzoic Acid (Calorimetric Standard)	Thianthrene Combustion Calorimetry	Synthetic Refuse-Derived Fuel	Nicotinic Acid (Combustion Calorimetric Standard)	Urea (Combustion Calorimetric Standard)
Unit of Issue	(30 g)	(30 g)	(100 g)	(25 g)	(25 g)
Heat of Combustion (in MJ/kg)*	26.434	33.480	13.87**	22.184	10.536

^{*} The calorific values (MJ/kg) may decrease upon the aging or normal oxidation of the coals. NIST will continue to monitor these calorific values and report any substantive change to the purchaser.

** Gross calorific value or HHV (Higher Heating Value).