

201.2 - Biological Buffer Systems (powder form)

These SRMs are used to calibrate clinical instruments (e.g., blood pH measurements), in the physiologically important range of pH 7 to 8. They are based on a biological buffer system for clinical pH measurements and are each certified for use as an admixture only. The pH(S) values for the buffer solutions are certified at 0.05 molal and 0.08 molal with respect to the free acid and the sodium salt admixture as a function of temperature. The certified temperature range is from 0° C to 50° C.

Technical Contact: kenneth.pratt@nist.gov

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM	2181	2182	2183	2184
Description	HEPES Free Acid	NaHEPESate	MOPSO Free Acid	NaMOPSOate
Unit Size	(60 g)	(60 g)	(50 g)	(50 g)
pH(S) Values (at 37 C)				
0.05 molal	7.364		6.699	
0.08 molal	7.373		6.694	