

105.2 - Serum Materials (frozen, liquid, and lyophilized forms)

These SRMs serve a variety of clinical measurement needs. SRM 909b is a lyophilized human serum for use in determining specified constituents. SRM 927d is a bovine serum albumin in a sterile 7% solution for use in the calibration and standardization of procedures to analyze total serum protein. SRM 956c is a frozen human serum for use in the calibration and standardization of procedures for the determination of specific electrolytes in either diluted or undiluted human serum or plasma. SRM 965b is a frozen human serum for evaluating the accuracy of procedures used to determine glucose in human serum and to validate secondary reference materials. SRM 967a is a frozen human serum for evaluating the accuracy of procedures for the determination of creatinine in human serum. SRM 968d is a lyophilized human serum for validating methods used to determine fat-soluble vitamins, carotenoids, and cholesterol in human serum and plasma. SRM 970 is a frozen human serum for validating methods for determining ascorbic acid in human serum and similar matrices. SRM 971 is a hormone in frozen human serum for evaluating the accuracy of procedures for the determination of the steroid hormones cortisol and progesterone in human serum. SRM 972 is a frozen human serum for evaluating the accuracy of clinical procedures for the determination of total cholesterol, HDL-cholesterol, LDL-cholesterol and triglycerides (triglycerides and total glyceride species). SRM 995b is a frozen human serum for evaluating vitamin D metabolites. SRM 995b is a frozen human serum for evaluating the accuracy of clinical procedures for the determination of total cholesterol. SRM 995b is a freeze-dried human serum for evaluating PCB congeners, chlorinated pesticides and total cholesterol in human serum and similar matrices. SRMs 995b and 995b are freeze-dried human serums for evaluating PCB congeners, chlorinated pesticides, and PBDE congeners with non-certified values for PFCs and dioxins/furans. SRM 2921 is a human cardiac troponin complex. SRM 2921 is primarily intended for use in calibrating clinical procedures and devices for the determination of cardiac troponin I (cTnI) in human serum. It can also be used for value-assignment of calibrators and control materials.

For further information see: [SP 260.83](#)

Technical Contact: karen.phinney@nist.gov

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM	Description	Unit Size	Analyte	909b-1	909b-2	BSA Concentration	Level I	Level II	Level III	Level IV	Single Concentration
909b	Human Serum	(2x3 each conc)+(6 water)	Calcium Chloride Cholesterol Creatinine Lithium Magnesium Potassium Sodium Total Glycerides Glycerides Urea	2.218 mmol/L 89.11 mmol/L 3.787 mmol/L 0.05618 mmol/L 0.6145 mmol/L 0.7634 mmol/L 3.424 mmol/L 120.76 mmol/L 0.949 mmol/L 0.804 mmol/L 5.51 mmol/L	3.532 mmol/L 119.43 mmol/L 6.084 mmol/L 0.4674 mmol/L 2.600 mmol/L 1.918 mmol/L 6.278 mmol/L 141.0 mmol/L 1.529 mmol/L 1.271 mmol/L 30.75 mmol/L						
927d	Bovine Serum Albumin (7%, solution)	set (10)				65.41 mg/L					
956c	Electrolytes in Frozen Human Serum	6 ampoules x 2.0mL each	CalciumIonized Calcium Chloride Lithium Magnesium Potassium Sodium				2.981 mmol/L 1.78 mmol/L 104.9 mmol/L 1.606 mmol/L 1.247 mmol/L 5.976 mmol/L 118.8 mmol/L	2.538 mmol/L 1.48 mmol/L 121.5 mmol/L 1.068 mmol/L 0.857 mmol/L 3.977 mmol/L 137.5 mmol/L	2.095 mmol/L 1.19 mmol/L 137.4 mmol/L 0.457 mmol/L 0.470 mmol/L 1.982 mmol/L 157.4 mmol/L		
965b	Glucose in Frozen Human Serum	set (8) (2 each conc)					1.836 mmol/L	4.194 mmol/L	6.575 mmol/L	16.35 mmol/L	
967a	Creatinine in Frozen Human Serum	set(4) (2 each conc)	Creatinine				0.0749 mmol/L	0.3427 mmol/L			
968d	Fat-Soluble Vitamins,Carotenoids and Cholesterol in Human Serum	set (2) (single conc)	Total-Retinol gamma-Tocopherol α-Tocopherol Total β-Carotene Cholesterol				8.41 μmol/L	28.05 μmol/L			1.09 μmol/L 3.34 μmol/L 13.77 μmol/L 0.145 μmol/L 3453 μmol/L
970	Ascorbic Acid in Frozen Human Serum	set (4) (2 each conc)	Ascorbic Acid				10.07 μmol/L	30.52 μmol/mol			
971	Hormones in Frozen Human Serum	5 mL	Male: 1 bottle Cortisol 297.0 nmol/L Progesterone 0.129 nmol/L Female: 1 bottle Cortisol 249.5 nmol/L Progesterone 6.11 nmol/L								
972	Vitamin D in Human Serum	set (4) (1 each level)	25-hydroxvitamin D2 25-hydroxvitamin D3 3-epi-25hydroxyvitamin D3				(1.46 nmol/L) 59.6 nmol/L (3.46 nmol/L)	4.14 nmol/L 30.8 nmol/L (1.9 nmol/L)	64.1 nmol/L 46.2 nmol/L (2.65 nmol/L)	5.81 nmol/L 82.3 nmol/L 94.1 nmol/L	
1951b	Lipids in Frozen Human Serum	set (4) (2 each conc)	Total Cholesterol Total Glycerides Triglycerides Only				4.804 mmol/L 1.370 mmol/L 1.208 mmol/L	6.895 mmol/L 2.988 mmol/L 2.700 mmol/L			
1955	Homocysteine and Folate in Frozen Human Serum	set (3) (1 each conc)	Homocysteine (μmol/L) 5-Methyltetrahydrofolic Acid Folic Acid (nmol)				3.98 μmol/L 4.26 nmol/L (0.49 nmol/L)	8.85 μmol/L 9.93 nmol/L (1.05 nmol/L)	17.7 μmol/L 37.1 nmol/L (1.07 nmol/L)		
1957	Organic Contaminants in Non-Fortified Human Serum	5 vials x 10 mL	PCB Congeners Chlorinated Pesticides PBDE Congeners								
1958	Organic Contaminants in Fortified Human Serum	5 vials x 10 mL	PCB Congeners Chlorinated Pesticides PBDE Congeners PBB 153 Congener (PFCs) (Dioxins/Furans)								
2921	Human Cardiac Troponin Complex	5x115 uL	cTnl cTnT cTnC				31.2 mg/L (36.9 mg/L) (24.2 mg/L)				

Values in parenthesis are not certified and given as reference values

105.2 - Serum Materials (frozen, liquid, and lyophilized forms)

These SRMs serve a variety of clinical measurement needs. SRM [900b](#) is a lyophilized human serum for use in determining specified constituents. SRM [907a](#) is a bovine serum albumin in a sterile 7% solution for use in the calibration and standardization of procedures to analyze total serum protein. SRM [956a](#) is a frozen human serum for use in the calibration and standardization of procedures for the determination of specific electrolytes in either diluted or undiluted human serum or plasma. SRM [957a](#) is a frozen human serum for evaluating the accuracy of procedures used to determine glucose in human serum and to validate secondary reference materials. SRM [962a](#) is a frozen human serum for evaluating the accuracy of procedures for the determination of creatinine in human serum. SRM [968a](#) is a lyophilized human serum for validating methods used to determine fat-soluble vitamins, carotenoids, and cholesterol in human serum and plasma. SRM [970a](#) is a frozen human serum for validating methods for determining ascorbic acid in human serum and similar matrices. SRM [971a](#) is a hormone in frozen human serum for evaluating the accuracy of procedures for the determination of the steroid hormones cortisol and progesterone in human serum. SRM [974a](#) is a frozen human serum for evaluating the accuracy of clinical procedures for the determination of total cholesterol, HDL-cholesterol, LDL-cholesterol and triglycerides (triglycerides and total glyceride species). SRM [977a](#) is a frozen human serum for evaluating vitamin D metabolites. SRM [981a](#) is a frozen human serum for evaluating the accuracy of clinical procedures for the determination of total cholesterol. SRM [988a](#) is a freeze-dried human serum for evaluating PCB congeners, chlorinated pesticides and total cholesterol in human serum and similar matrices. SRMs [9857](#) and [9858](#) are freeze-dried human serums for evaluating PCB congeners, chlorinated pesticides, and PBDE congeners with non-certified values for PFCs and dioxins/furans. SRM [9921](#) is a human cardiac troponin complex. SRM 2921 is primarily intended for use in calibrating clinical procedures and devices for the determination of cardiac troponin I (cTnI) in human serum. It can also be used for value-assignment of calibrators and control materials.

For further information see: [SP 260.83](#)
Technical Contact: karen.phinney@nist.gov

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

Values in parenthesis are not certified and given as reference values