105.6 - Drugs of Abuse in Urine, Multianalyte (powder form)

These SRMs are for verifying the accuracy of methods used to determine cocaine, morphine and codeine and opiate-based substances classed as drugs of abuse in humans. Each SRM consists of multianalyte drugs of abuse. SRM 1511 consists of a mixture of five substances -- morphine, codeine, cocaine metabolite, marijuana metabolite and phencyclidine in freeze-dried urine and is provided as a set of three bottles, each containing all analytes (no blank). SRM 2381 consists of freeze-dried urine and is provided in a set of four 20 mL bottles -- three levels plus a blank.

<u>Table 105.5</u> lists the Single Analyte Drugs of Abuse in Urine.

Technical Contact: <u>karen.phinney@nist.gov</u>

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

| SRM | Description | Unit of Issue | Benzoylecgonine | Morphine | Codeine | THC-9-COOH | Phencyclidine |
|------|-------------------------------|---------------|-----------------|----------|---------|------------|---------------|
| 1511 | Multi Drugs of Abuse in Urine | 3 bottles | 162 | 309 | 288 | 14.1 | 23.8 |
| 2381 | Morphine and Codeine in Urine | set (4) | | 134-580 | 130-560 | | |

105.6 - Drugs of Abuse in Urine, Multianalyte (powder form)

These SRMs are for verifying the accuracy of methods used to determine cocaine, morphine and codeine and opiate-based substances classed as drugs of abuse in humans. Each SRM consists of multianalyte drugs of abuse. SRM 1511 consists of a mixture of five substances -- morphine, codeine, cocaine metabolite, marijuana metabolite and phencyclidine in freeze-dried urine and is provided as a set of three bottles, each containing all analytes (no blank). SRM 2381 consists of freeze-dried urine and is provided in a set of four 20 mL bottles -- three levels plus a blank.

<u>Table 105.5</u> lists the Single Analyte Drugs of Abuse in Urine.

Technical Contact: <u>karen.phinney@nist.gov</u>

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.